



# 1. Contraindications

## 1.1. General

Do not treat patients if the following conditions are present:

- Bacterial or viral infection, acute inflammations
- Impaired immune system
- Isotretinoin in the past 12 months
- Skin-related autoimmune diseases
- Radiation therapy and chemotherapy
- Poor healing and unhealed wounds in the treatment area
- Metal implants near treated area or neutral electrode
- Permanent implant near the area
- Pacemaker or internal defibrillator, or any other active electrical implant anywhere in the body
- Current condition or history of skin cancer, or current condition of any other type of cancer, or pre-malignant moles
- History of any type of malignant cancer
- Active collagen diseases
- Cardiovascular diseases (such as vascular diseases, peripheral arterial disease, thrombophlebitis and thrombosis)
- Pregnancy/nursing or IVF procedure
- History of bleeding coagulopathies, use of anticoagulants
- Any active condition in the treatment area, such as sores, psoriasis, eczema, rash and rosacea
- Any surgical procedure in the treatment area within the last three months or before complete healing
- Poorly controlled endocrine disorders, such as diabetes
- Tuberculosis
- Hepatitis
- Febrile conditions
- Acute neuralgia and neuropathy
- Kidney or liver failure
- Sensitivity disorders in the treatment area
- Varicose veins, pronounced edemas
- Skin dermabrasion, skin resurfacing, or deep chemical peeling in the treatment area within 3 months prior to the treatment

For the full range of contraindications, warnings, and cautions, refer to the Operator's Manual.